# **OMRON**

# GETTING **STARTED**

**Premium Arm Blood Pressure** Monitor

BP755

with ComFit<sup>™</sup> Cuff



The American Heart Association's guideline for upper limit of normal home blood pressure is:


**135** <sup>mmHg</sup> systolic **85** <sup>mmHg</sup> diastolic

Please read the instruction manual for complete product information.

#### Have questions? **DO NOT RETURN THIS PRODUCT TO THE STORE.**

Call our Toll-Free Consumer Help Line at 800-634-4350 for assistance with your Omron product.

www.omronhealthcare.com









- To Bypass Date and Time Setup Go To **B**.
- **B.** Press **START/STOP** button to turn unit on. The year will flash. Press **START/STOP** again to turn unit off. You have now bypassed setting the Date and Time. Go to **STEP 3.**

For more information on date and time setting, see pages 17-19 in the instruction manual.


A. Press **START/STOP** button once. The year will flash. Press the Memory Up button 🜔 to advance to current year, then press 🦲 to confirm. (Press the Memory Down button to go back to the previous value.) Change month, day and time the same way. Unit will shut off when complete.





Place cuff directly on bare skin of upper left arm, with cuff edge 1/2 inch above elbow.

Align cuff so the blue stripe is centered on the inside of your arm with the tubing pointing down toward your hand. Close cuff snugly.



To use on right arm, see page 22 in instruction manual.

### STEP To turn unit on and take a measurement, press **START/STOP** button once.

The cuff will inflate automatically. Hold still until cuff deflates and results are displayed.



### **FREQUENTLY ASKED QUESTIONS**

## **Q.** Why are my blood pressure readings sometimes different?

**A.** Your blood pressure changes constantly based on your daily activities and time of day. Because of this, some fluctuations can be expected. That's why it is important to use averages over time to compare readings. Because each person is different, it's important to consult with your doctor about what your numbers mean to you.

#### **Q.** Why could my doctor's reading be different than my readings taken at home?

**A.** Having your blood pressure measured by a healthcare professional in a doctor's office can cause nervousness and may result in an elevated reading. Because these readings can vary considerably, you should take your blood pressure on a regular basis at home and share your results with your healthcare provider.

### **OTHER IMPORTANT TIPS:**

- For accurate readings, sit quietly with feet flat on the floor and rest for five minutes before taking measurements. Do not measure within 30 minutes after eating, exercising, bathing, smoking or drinking alcohol.
- To review readings, press the 🔾 and 🜔 memory buttons to scroll through readings.
- To use the Advanced Averaging function and display averages, see page 26 in the instruction manual.
- To keep your monitor in the best condition, follow the Care and Maintenance directions on page 29 of the instruction manual.

