
OMRON®

INSTRUCTION MANUAL

AUTOMATIC INFLATION BLOOD PRESSURE MONITOR

Model **HEM-712CLC**

With Large Adult Cuff

ESPAÑOL

ENGLISH

TABLE OF CONTENTS

Introduction	3
Important Safety Notes	4
Know Your Unit	5
Suggestions Before Taking A Measurement	6
Battery Installation/Replacement	7
Care and Maintenance	8
How To Apply The Arm Cuff	9-10
How To Take A Measurement	11-12
How To Use The Memory Function	13
How To Use The AC Adapter	14
Error Indicators	15
FCC Statement	16
Warranty Information	17-18
Specifications	19

INTRODUCTION

Thank you for purchasing the Omron HEM-712CLC Automatic Inflation Blood Pressure Monitor.

Fill-in for future reference

DATE PURCHASED: _____

SERIAL NUMBER: _____

Staple your purchase receipt here.

This monitor features automatic inflation and deflation to quickly measure your blood pressure and pulse rate. The measurement reading is clearly displayed on a large digital panel.

The HEM-712CLC digital blood pressure monitor uses the oscillometric method of blood pressure measurement. This means the monitor detects your blood's movement through your brachial artery and converts the movements into a digital reading. An oscillometric monitor does not need a stethoscope so the monitor is simple to use.

The Omron HEM-712CLC Automatic Inflation Blood Pressure Monitor is intended for home use.


SAVE THESE INSTRUCTIONS

IMPORTANT SAFETY NOTES

To assure the correct use of the product basic safety measures should always be followed including the precautions listed below:

- Read all information in the instruction book and any other literature included in the box before using the unit.
- Contact your physician for specific information about your blood pressure. Follow the instructions of your healthcare provider.
- Operate the unit only as intended. Do not use for any other purpose.
- The unit is intended for use in measuring blood pressure and pulse rate in adult patient population. Do not use on infants.
- Do not use a cellular phone near the unit. It may result in operational failure.
- Changes or modifications not approved by Omron Healthcare will void the user warranty. Do not disassemble or attempt to repair the unit or components.

KNOW YOUR UNIT


SUGGESTIONS BEFORE TAKING A MEASUREMENT


- 1.** Avoid eating, smoking, and exercising for 30 minutes before taking a measurement. Rest for at least 15 minutes before taking the measurement.
- 2.** Stress raises blood pressure. Avoid taking measurements during stressful times.
- 3.** Remove tight-fitting clothing from your left arm.
- 4.** Measurements should be taken in a quiet place.
- 5.** Sit in a chair with your feet flat on the floor. Rest your left arm on a table so that the cuff is at the same level as your heart.
- 6.** Remain still and do not talk during the measurement.
- 7.** Keep a record of your blood pressure and pulse readings for your physician. A single measurement does not provide an accurate indication of your true blood pressure. You need to take and record several readings over a period of time. Try to measure your blood pressure at the same time each day for consistency.
- 8.** Wait 2-3 minutes between measurements. The wait time allows the arteries to return to the condition prior to taking the blood pressure measurement. You may need to increase the wait time depending on your individual physical characteristics.

BATTERY INSTALLATION / REPLACEMENT


1. Slide the battery cover off in the direction of the arrow.


2. Install 4 "AA" size batteries so + (positive) and - (negative) polarities match the polarities of the battery compartment as indicated.

3. Replace the battery cover.

NOTE: If the Low Battery Indicator appears on the display, replace all four batteries. Long-life alkaline batteries are recommended.

CARE AND MAINTENANCE

To keep your digital blood pressure monitor in the best condition and protect the unit from damage follow the directions listed below:

Clean the monitor with a soft dry cloth. Do not use any abrasive or volatile cleaners. Do not attempt to clean the cuff. Never immerse the monitor or any components in water.

Store the monitor in a safe and dry location. Do not fold the cuff and tubing tightly. Do not subject the monitor to extreme hot or cold temperatures, humidity and direct sunlight.

Avoid subjecting the monitor to strong shocks, such as dropping the unit on the floor.


Remove the batteries if the unit will not be used for three months or longer. Always replace all the batteries with new ones at the same time.

Use the unit consistent with the instruction provided in this manual. Use only authorized parts and accessories.

HOW TO APPLY THE ARM CUFF


1. If the cuff is assembled correctly, the sewn hook material will be on the outside of the cuff loop and the metal D-ring will not touch your skin.

2. Pass the end of the cuff furthest from the tubing through the metal D-ring to form a loop. The smooth cloth should be on the inside of the cuff loop.


3. Put your left arm through the cuff loop. The bottom of the cuff should be approximately 1/2 inch above the elbow. The green marker on the cuff should lie over the brachial artery on the inside of the arm. The cuff tube should run down center of arm even with middle finger.

1/2 inch (12.7 mm)




DO NOT BLOCK THE
AIR RELEASE VALVE

NOTE: Do not block the end of the short tube which contains a supplemental air release device that is necessary to maintain a proper deflation ratio. Blocking this tube may result in an inaccurate reading.

HOW TO APPLY THE ARM CUFF

4. Pull the cuff so that the top and bottom edges are tightened evenly around your arm.
5. When the cuff is positioned correctly, press the sewn hook material FIRMLY against the pile side of the cuff.
6. Make certain the cuff fits snugly around your arm. The cuff should make good contact with your skin.
7. Sit in a chair with your feet flat on the floor and place your arm on a table so that the cuff is at the same level as your heart.
8. Relax your arm and turn your palm upward.
9. Be sure there are no kinks in the air tubing.


HOW TO TAKE A MEASUREMENT

Be sure to read “How to Apply the Arm Cuff” before taking a reading.

1. Press the ON/OFF button.

- a) All display symbols appear for approximately one second. This is the initial LCD test.
- b) The display symbols disappear and the Deflation Indication symbol (▼) starts to flash.
- c) When the monitor completes the necessary preparations before measurement, the Heart Symbol (♥) appears in the screen.

Initial
LCD Test

2. By pressing the START button the unit automatically inflates the cuff to 170 mmHg.

- a) If during your measurement the monitor determines that a higher pressure value is required, it will automatically increase the pressure value by 40 mmHg.
- b) To inflate higher than 210 mmHg, press and hold the START button until the desired pressure value is reached (approximately 40 mmHg higher than your estimated systolic blood pressure).
Do not inflate higher than 280 mmHg.

Deflation
Indication


Ready
to Measure


NOTE: If you want to stop the measurement, push the ON/OFF button. The monitor will stop inflating and start deflating rapidly, then the monitor will power off.

HOW TO TAKE A MEASUREMENT

3. When the correct pressure value is reached, the cuff starts to deflate automatically. As the cuff deflates, decreasing numbers appear on the display. The Heart Symbol (♥) flashes at every heartbeat.

4. When the measurement is complete, the Heart Symbol (♥) stops flashing and your blood pressure and pulse rate are displayed alternately. Results will be displayed until you turn the unit off. The unit will automatically shut off in approximately 5 minutes.


NOTE: Wait 2-3 minutes before taking another blood pressure measurement. Waiting between readings allows the arteries to return to the condition prior to taking the blood pressure measurement.

HOW TO USE THE MEMORY FUNCTION


This monitor has a memory capable of storing fourteen (14) sets of readings. Every time you complete the measurement, the monitor automatically stores the blood pressure reading and pulse rate.

- Push the MEMORY button to recall stored readings while the Heart Symbol (♥) is displayed.


- When fourteen (14) sets of readings are stored in memory, the oldest set will be deleted to store a new set.
- To delete all stored readings, press and hold the MEMORY button, then press START button for 2 seconds. **All the readings will be deleted.**

NOTE: You cannot delete the stored readings partially.

HOW TO USE THE AC ADAPTER

Use only the authorized Omron AC Adapter, Model Number HEM-ADPT1, with this monitor. Use of any other adapter may cause damage to the monitor and will void the user warranty.


To connect the AC Adapter

- ① Insert the AC Adapter Plug into the AC Adapter Jack on the right side of the monitor as illustrated.
- ② Plug the AC Adapter into a 120V AC outlet (50-60 cycles).

To disconnect the AC Adapter remove the AC Adapter from the electrical outlet first and then remove the AC Adapter Plug from the monitor.

NOTE: The monitor is designed not to draw power from the batteries when the AC Adapter is used.

ERROR INDICATORS

Error Symbol	Cause	Correction
	Cuff under-inflated.	Remove arm cuff. Wait 2-3 minutes and take another measurement. Repeat steps listed under "How To Take A Measurement". Inflate the cuff 40 mmHg higher than the previous inflation value.
	Movement during Measurement. No pulse or blood pressure reading due to very weak artery pulse.	
	Cuff over-inflated (more than 300 mmHg).	Carefully read and repeat steps listed under "How To Take A Measurement."
	Batteries have run down.	Replace all four "AA" batteries.

FCC STATEMENT

Note:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)
This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE

(for Canada only)

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the interference-causing equipment standard entitled " Digital Apparatus", ICES-003 of the Canadian Department of Communications.

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Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

WARRANTY INFORMATION

Limited Warranties

Your HEM-712CLC Automatic Blood Pressure Monitor, excluding the monitor cuff, is warranted to be free from defects in materials and workmanship appearing within 2 years, when used in accordance with the instructions provided with the monitor. The monitor cuff is warranted to be free from defects in materials and workmanship appearing within one year, when the monitor is used in accordance with the instructions provided with the monitor. The above warranties extend only to the original retail purchaser.

We will, at our option, repair or replace without charge any monitor or monitor cuff covered by the above warranties. Repair or replacement is our only responsibility and your only remedy under the above warranties.

To obtain warranty service, ship the monitor and all of the components, together with proof of purchase and a note explaining the problem with \$5.00 for return shipping and insurance to:

Omron Healthcare, Inc.
Attn: Repair Department
300 Lakeview Parkway
Vernon Hills, IL 60061

Be sure to include the model number, your name and your phone number on the enclosed note and any other correspondence.

ALL IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE APPLICABLE WRITTEN WARRANTY ABOVE. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

OMRON SHALL NOT BE LIABLE FOR LOSS OF USE OR ANY OTHER INCIDENTAL, CONSEQUENTIAL OR INDIRECT COSTS, EXPENSES OR DAMAGES. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusions may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

FOR CUSTOMER SERVICE

Visit our web site at:

www.omronhealthcare.com

Call toll free:

1-800-634-4350

WARRANTY INFORMATION

Limited Warranties (Outside of the USA and Canada)

Your HEM-712CLC Automatic Blood Pressure Monitor, excluding the monitor cuff, is warranted to be free from defects in materials and workmanship appearing within 2 years, when used in accordance with the instructions provided with the monitor. The monitor cuff is warranted to be free from defects in materials and workmanship appearing within one year, when the monitor is used in accordance with the instructions provided with the monitor. The above warranties extend only to the original retail purchaser.

Should repair be needed during the warranty period, please contact the authorized distributor in your country.

The above warranty is complete and exclusive. The warrantor expressly disclaims liability for incidental, special, or consequential damages of any nature. (Some countries do not allow the exclusion or limitation of incidental or consequential damages, so the above warranty may not apply to you.)

Any implied warranties arising by the operation of law shall be limited in duration to the term of this warranty. (Some countries do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.)

This warranty gives you specific legal rights and you may have other rights which vary country to country.

SPECIFICATIONS

Model:	HEM-712CLC
Display:	LCD Digital Display
Measurement Range:	Pressure: 0 to 280 mmHg Pulse: 40-180/min.
Accuracy/Calibration:	Pressure: ± 3 mmHg or 2% of reading Pulse: $\pm 5\%$ of reading
Inflation:	Automatic by electric pump
Deflation:	Automatic pressure release valve
Rapid Air Release:	Automatic exhaust valve
Pressure Detection:	Capacitive pressure sensor
Measurement Method:	Oscillometric method
Pulse Wave Detection:	Capacitance type pressure sensor
Power Source:	4 "AA" batteries (not included), AC Adapter (model HEM-ADPT1)
Battery Life:	Approx. 300 uses when used once a day for 2 minutes (Usage is decreased when using the large cuff.)
Operating Temperature/ Humidity:	50°F to 104°F (10°C to 40°C) 30 to 85% RH maximum
Storage Temperature/ Humidity:	-4°F to 140°F (-20°C to 60°C) 10 to 95% RH maximum
Console Weight:	Approximately 13 oz. (350g) with cuff, excluding batteries
Outer Dimensions:	Approx. 5" (l) x 4 3/4" (w) x 3 3/4" (h) (125mm x 119mm x 94mm)
Cuff Dimensions:	Approx. 6 7/8" (w) x 24 1/2" (l) (175mm x 622mm)
Arm Circumference Range:	Approximately 13" to 17" (330mm - 432mm)
Contents:	Large Arm Cuff, Monitor, Instruction Manual, and AC Adapter
UPC Code:	0 73796 07124 0

NOTE: These specifications, to improve performance, are subject to change without notice.

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Copyright © 2003 Omron Healthcare, Inc.
© 2003 Omron Healthcare, Inc.