



Body Composition Monitor and Scale Comparison Chart



Model Number	Body Composition Monitor and Scale HBF-514C	Body Composition Monitor and Scale BCM-500
Body Weight Range	0 to 330 lbs. with an increment of .2 lbs.	4.4 to 330 lbs. with an increment of .2 lbs.
Body Fat Percentage	5.0 to 60.0% with an increment of .1%	5.0 to 60.0% with an increment of .1%
Body Mass Index	7.0 to 90.0 with an increment of 0.1	7.0 to 90.0 with an increment of 0.1
Skeletal Muscle Percentage	5.0 to 50.0% with an increment of .1%	5.0 to 50.0% with an increment of .1%
Visceral Fat Level	30 levels with an increment of 1 level	30 levels with an increment of 1 level
Resting Metabolism	385 to 3999 kcal with an increment of 1 kcal	385 to 3999 kcal with an increment of 1 kcal
Body Age	18 to 80 year old with 1 year increments	N/A
Memory	90-Day	120 Total Readings in Unit, Unlimited in OMRON connect US/CAN app
Bluetooth		✓
User Profiles	4	4
Sensor Technology	Full Body	Foot - to - Foot
Battery	4 AA	4 AAA