## **OMRON Body Composition Monitor and Scale Comparison Chart** 1382 **Body Composition Monitor and Scale Body Composition Monitor and Scale Digital Scale with Product Name** with Bluetooth® Connectivity with Seven Fitness Indicators **Bluetooth® Connectivity** and Model Number **HBF-514C BCM-500** SC-150 **Body Weight Range** 0 to 330 lbs. with an increment of .2 lb. 4.4 to 330 lbs. with an increment of .2 lb. 11 to 330 lbs. with an increment of .2 lb **Body Fat Percentage** 5.0 to 60.0% with an increment of .1% 5.0 to 60.0% with an increment of .1% N/A **Body Mass Index** 7.0 to 90.0 with an increment of 0.1 7.0 to 90.0 with an increment of 0.1 N/A **Skeletal Muscle Percentage** 5.0 to 50.0% with an increment of .1% 5.0 to 50.0% with an increment of .1% N/A Visceral Fat Level 30 levels with an increment of 1 level 30 levels with an increment of 1 level N/A **Resting Metabolism** N/A 385 to 3999 kcal with an increment of 1 kcal 385 to 3999 kcal with an increment of 1 kcal **Body Age** 18 to 80 year old with 1 year increments N/A N/A 30 Readings for 4 Users 30 Total Readings for One User Memory 90-Day (Up to 120 Total Readings) $\checkmark$ **Bluetooth** N/A 4 **User Profiles** 1 **Sensor Technology** Full Body Foot - to - Foot N/A 4 AAA 4 AAA **Battery** 4 AA