OMRON.	Upper Arm Blood Pressure Monitor Comparison Chart										
Model Number	3 Series® BP7100	Bronze BP5100	5 Series® BP7200	5 Series® Wireless BP7250	Silver BP5250	7 Series® Wireless BP7350	Gold BP5350	10 Series [®] Wireless BP7450	Platinum BP5450	Evolv® BP7000	Complete™ BP7900
Advanced Accuracy	√	√	√	√	√	√	√	√	√	√	√
Hypertension Indicator			√	\checkmark	\checkmark	~	~	√	√		
Bluetooth®				\checkmark	\checkmark	~	~	√	√	√	√
TruRead								1	√		
Extra-Large Digits w/Backlight								1	√		
One-Touch Use	\checkmark									√	
Guest Mode				√		~	\checkmark	√	√		
Two User Mode			√			~	\checkmark	√	√		
Body Movement Detection	\checkmark	1	√	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	√	~	
Advanced Averaging			√	1	1	~	~	√	√		
Irregular Heartbeat Detector	\checkmark	1	√	√	√	~	~	1	√	√	√
Device Storage Capacity Memory per User	14	14	60	60	80	60 per user	60 per user	100 per user	100 per user	100	90
Morning Hypertension Indicator							~		√		
Morning Average/Week							~		√		
EKG											√
Side-by-Side Comparison Display						√	√	√	√		
Cuff	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"
Replacement Cuff	CD-WR17	CD-WR17	HEM-RML31-B	HEM-RML31-B	HEM-FL31-B	HEM-FL31-B	HEM-RML31-B	HEM-FL31-B	HEM-RML31-B		HEM-RML31-B
Optional Small Cuff (7"- 9")	CD-CS9	CD-CS9	HEM-CS24-B	HEM-CS24-B			HEM-CS24-B		HEM-CS24-B		HEM-CS24-B
Storage Case					√		\checkmark		\checkmark		\checkmark
Power Options	4 AA Not Included / Not compatible with AC Adaptor	4 AA Not Included / Not compatible with AC Adaptor	4 AA Included / AC Adapter Not Included	4 AA Included / AC Adapter Included	4 AA Included / AC Adapter Included	4 AAA Included / Not compatible with AC Adaptor	4 AA Included / Not compatible with AC Adaptor				

May-23