














### Upper Arm Blood Pressure Monitor Comparison Chart

											
Model Number	3 Series® BP7100	Bronze BP5100	5 Series® BP7200	5 Series® Wireless BP7250	Silver BP5250	7 Series® Wireless BP7350	Gold BP5350	10 Series® Wireless BP7450	Platinum BP5450	Evolv® BP7000	Complete™ BP7900
Advanced Accuracy	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Hypertension Indicator			✓	✓	✓	✓	✓	✓	✓		
Bluetooth®				✓	✓	✓	✓	✓	✓	✓	✓
TruRead								✓	✓		
Extra-Large Digits w/Backlight								✓	✓		
One-Touch Use	✓									✓	
Guest Mode			✓	✓		✓	✓	✓	✓		
Two User Mode			✓			✓	✓	✓	✓	✓	
Body Movement Detection	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Advanced Averaging			✓	✓	✓	✓	✓	✓	✓		
Irregular Heartbeat Detector	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Device Storage Capacity Memory per Morning Hypertension Indicator	14	14	60	60	80	120 (60 per user)	120 (60 per user)	200 (100 per user)	200 (100 per user)	100	90
Morning Average/Week							✓		✓		
EKG											✓
Side-by-Side Comparison Display						✓	✓	✓	✓		
Cuff	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"	Easy-Wrap ComFit 9"-17"	Wide-Range D-Ring 9"-17"
Replacement Cuff	CD-WR17	CD-WR17	HEM-RML31-B	HEM-RML31-B	HEM-FL31-B	HEM-FL31-B	HEM-RML31-B	HEM-FL31-B	HEM-RML31-B		HEM-RML31-B
Optional Small Cuff (7"-9")	CD-CS9	CD-CS9	HEM-CS24-B	HEM-CS24-B			HEM-CS24-B		HEM-CS24-B		HEM-CS24-B
Storage Case					✓		✓		✓		✓
Power Options	4 AA Not Included / Not compatible with AC Adaptor	4 AA Not Included / Not compatible with AC Adaptor	4 AA Included / AC Adaptor Not Included	4 AA Included / AC Adaptor Not Included	4 AA Included / AC Adaptor Not Included	4 AA Included / AC Adaptor Not Included	4 AA Included / AC Adaptor Not Included	4 AA Included / AC Adaptor Included	4 AA Included / AC Adaptor Included	4 AAA Included / Not compatible with AC Adaptor	4 AA Included / Not compatible with AC Adaptor