	Wrist Blood Pressure Monitor Comparison		
OMRON			
Model Number	3 Series® BP6100	7 Series® BP6350	Gold BP4350
UltraSilent Inflation		√	\checkmark
Cuff	5.3"-8.5"	5.3"-8.5"	5.3"-8.5"
Hypertension Indicator	\checkmark	\checkmark	\checkmark
Bluetooth ®		\checkmark	\checkmark
Two User Mode			\checkmark
Advanced Averaging	\checkmark	\checkmark	\checkmark
Irregular Heartbeat Detector	\checkmark	\checkmark	\checkmark
Total Memory	60	90	200
Storage Case	\checkmark	\checkmark	\checkmark
Slim & Portable Design	\checkmark	\checkmark	\checkmark
Heart Zone Guidence		\checkmark	\checkmark
Morning Hypertension Indicator			\checkmark
Morning Average/Week			\checkmark
Power Options	2 AAA Batteries Included	2 AAA Batteries Included	2 AAA Batteries Included