

Complete™ Blood Pressure Monitor and EKG Device (BP7900)

Complete™ from OMRON is the **first blood pressure monitor with EKG capability in a single device for home use in the U.S.**, providing a more complete heart health picture.

- An at-home EKG and blood pressure monitoring device from the #1 doctor and pharmacist recommended blood pressure monitor brand, **Complete represents groundbreaking innovation for millions with AFib or a family history of irregular heartbeat.** Launched in the U.S. in 2019, Complete is registered by the FDA as a medical device.
- Up to six million Americans have atrial fibrillation (AFib) – a condition that increases their risk of stroke five-fold. Complete has the **unprecedented ability to monitor two critical risk factors for stroke** – AFib and high blood pressure – at the same time, in the patient’s home, and with a single device.
- **Complete syncs with the OMRON Connect mobile app** in order to provide Intuitive and user-friendly integration. Blood pressure and EKG readings can be stored, tracked, and shared with physicians in order to gain additional insight and improve treatment over time.
- For **easy access and quick measurement of EKG**, Complete features touchpad electrodes that are located on the top face and both sides of the monitor. Users can store up to 90 readings on the monitor and can transfer readings to the OMRON Connect app for unlimited storage in a secure, HIPAA-compliant cloud.
- With EKG an essential heart health measurement, **Complete uses an advanced algorithm designed by AliveCor to improve detection.** In tandem with trusted medical-grade blood pressure measurement from OMRON, the AliveCor technology alerts doctors to the possibility of AFib, bradycardia, tachycardia, and sinus rhythm.



About this Innovation

“We designed Complete to be a difference-maker for patients with AFib and to give anyone who suspects they may have irregular heartbeat a tool at home, available anytime, to check for this condition. We’ve heard from

numerous consumers and physicians that Complete is breaking through at achieving first detection of this dangerous condition,” said Carol Lucarelli, OMRON’s executive director of marketing and ecommerce.

Device Specifications

Monitor: Approximately 9.1" (w) × 3.9" (h) × 4.8" (l) (232 mm × 98 mm × 123mm)

Arm cuff: Approximately 5.7" × 23.4" (air tube: 29.5") (145 mm × 594 mm (air tube: 750 mm)

Availability & MSRP

OmronHealthcare.com, \$199.99; currently on sale for \$149.99 (USD)

Media Contact

MWW for OMRON Healthcare / Erin LaRossa | 201.316.3703 | elarossa@mww.com

More Information

<https://omronhealthcare.com/products/complete-upper-arm-blood-pressure-monitor-ekg-bp7900/>

**High-res images available on request*