



HeartGuide™ Wearable Blood Pressure Monitor (BP8000-M/BP8000-L)

HeartGuide is the **first wearable blood pressure monitor** in the compact form of a wristwatch – for those who want to keep close watch on their heart health anytime, anywhere and take action for their heart health.

- The sleek and stylish HeartGuide uses an inflatable oscillometric blood pressure cuff -- the same technology used in a doctor's office -- to ensure accurate and reliable readings.
- **HeartGuide received 510K FDA clearance as a medical device.**
- HeartGuide **accurately measures blood pressure, daily activity (including steps, distance and calories burned), pulse rate and sleep quality**. It also functions as a smartwatch, where users can receive text, call and email notifications.
- OMRON filed **more than 80 new patents** to create HeartGuide in the compact form of a wristwatch and provide medical-grade accuracy via wrist measurement. Some of the pumps and valves in HeartGuide are as small as a grain of rice.
- HeartGuide holds up to 100 measurements in memory and connects exclusively with the OMRON **HeartAdvisor™** app, which stores and tracks heart health data and provides insights and personalized coaching to help users take action for their heart health.
- A typical user can expect to charge HeartGuide™ approximately 2-3 times per week, depending upon the frequency of use of HeartGuide's feature.
- HeartGuide™ is available in sizes medium and large



Specifications

Diameter approximately 1.89" (48 mm)
Case thickness approximately 0.55" (14 mm)
Band width approximately 1.18" (30 mm)

Availability & MSRP

OmronHealthcare.com; \$499 (USD)



Media Contact

MWWPR for Omron Healthcare / Katherine Clark | 413.579.4791 | kclark@mww.com

More Information

<https://omronhealthcare.com/heartguide-wearable-blood-pressure-monitor-bp8000m/>

**High res images available on request*