



HeartGuide™ Wearable Blood Pressure Monitor (BP8000-M/L)

HeartGuide from OMRON® is the **first wearable blood pressure monitor**. It's made in the compact form of a wristwatch for those who want to keep close watch on their blood pressure anytime, anywhere, and take action for their heart health.

- HeartGuide uses an inflatable oscillometric blood pressure cuff -- the same technology used in a doctor's office -- to ensure accurate and reliable readings. Oscillometric technology is the FDA standard for accurate, automated, medical-grade personal blood pressure measurement.
- OMRON filed **more than 80 new patents** to create HeartGuide in the compact form of a wristwatch and provide medical-grade accuracy via wrist measurement. Some of the pumps and valves in HeartGuide are as small as a grain of rice.
- HeartGuide **accurately measures blood pressure, daily activity (including steps, distance and calories burned), pulse rate and sleep quality**. It also functions as a smartwatch, where users can receive text, call and email notifications.
- HeartGuide holds up to 100 measurements in memory and connects exclusively with the OMRON **HeartAdvisor™** app, which allows users to store unlimited data, track it over time, and share it with their doctor. HeartAdvisor provides insights and coaching to help users take real-time action for their heart health.
- A typical user can expect to charge HeartGuide approximately 2-3 times per week, depending upon the frequency of use.
- HeartGuide is now available in sizes medium and large in the U.S. at OmronHealthcare.com. HeartGuide recently launched in Japan and countries in the EU (including the U.K., Germany, France, Italy and Spain).



Specifications

Diameter approximately 1.89" (48 mm)
Case thickness approximately 0.55" (14 mm)
Band width approximately 1.18" (30 mm)



Availability & MSRP	OmronHealthcare.com; \$499 (USD)
--------------------------------	----------------------------------

More Information	https://omronhealthcare.com/heartguide-wearable-blood-pressure-monitor-bp8000m/
-------------------------	---

**High res images available on request*