

OMRON Nebulizer Comparison Chart

Find out which one is right for your child. Ask your healthcare provider for a prescription to complete your purchase.



**MicroAir®
Mesh Nebulizer**



**CompAir®
Compressor
Nebulizer**



**CompAir®
Compressor
Nebulizer for Kids**

NDC Number	73796-0451-00	73796-0458-01	73796-8001-47
Handheld	✓		
Tabletop		✓	✓
Kid-friendly design			✓
For children and adults	✓	✓	
Sound	<23 dBA; virtually silent operation	≈52dBA of sound while in use, similar to the sound of an electric toothbrush*	≈45dBA of sound while in use, a little more than a quiet library*
Power	AA Battery or AC Adapter	AC Adapter	AC Adapter
Weight	4.2 oz	9.5 oz; 3.5x lighter than most tabletop compressor models	6.3 oz; our lightest tabletop compressor model
360° operation for use in any position	✓		
Mass Mean Aerodynamic Diameter (MMAD)	5 μm	3 μm	3 μm
Warranty	2 years	3 years	3 years
Reimbursable	✓	✓	✓