

OMRON Nebulizer Comparison Chart

Find out which one is right for your child. Ask your healthcare provider for a prescription to complete your purchase.



**Portable MicroAir®
Battery-Operated
Nebulizer**



**CompAir®
Compressor
Tabletop
Nebulizer System**



**CompAir®
Compressor
Nebulizer with
Kid's Accessory**

| | | | |
|---|-------------------------------------|---|---|
| NDC Number | 73796-0451-00 | 73796-0458-01 | 73796-8001-47 |
| Handheld | ✓ | | |
| Tabletop | | ✓ | ✓ |
| Kid-friendly design | | | ✓ |
| For children and adults | ✓ | ✓ | |
| Sound | <23 dBA; virtually silent operation | ≈52dBA of sound while in use, similar to the sound of an electric toothbrush* | ≈45dBA of sound while in use, a little more than a quiet library* |
| Power | AA Battery or AC Adapter | AC Adapter | AC Adapter |
| Weight | 4.2 oz | 9.5 oz; 3.5x lighter than most tabletop compressor models | 6.3 oz; our lightest tabletop compressor model |
| 360° operation for use in any position | ✓ | | |
| Mass Mean Aerodynamic Diameter (MMAD) | 5 μm | 3 μm | 3 μm |
| Warranty | 2 years | 3 years | 3 years |
| Reimbursable | ✓ | ✓ | Pending |