

TENS App Compatible Devices & Operating Systems

Please note: You must have both the device and the operating systems below to be compatible with our mobile applications.

Apple		
Phone Model	Operating System	TENS
iPhone 5, 5c, 5S	iOS 9.0+	✓
iPhone 6, 6 Plus	iOS 9.0+	✓
iPhone 6S, 6S Plus	iOS 9.0+	✓
iPhone 7, 7Plus	iOS 10.0+	√
iPhone 8, 8+	iOS 11.0+	✓

Droid		
Phone Model	Operating System	TENS
Maxx 2	Android 6.0.1	✓
Turbo 2	Android 7.0	√

Google		
Phone Model	Operating System	TENS
Google Pixel	Android 8.0.0	√

HTC		
Phone Model	Operating System	TENS
HTC M9	Android 7.0	√
HTC A9	Android 6.0	√
HTC 10	Android 7.0	√

Huawei		
Phone Model	Operating System	TENS
Nexus 6P	Android 8.0.0	✓

updated 3/5/2021

LG		
Phone Model	Operating System	TENS
LG G5	Android 6.0.1	✓
LG G4	Android 6.0	✓
LG V20	Android 7.0	✓
Nexus 5	Android 6.0	✓
Nexus 5X	Android 7.1.2	✓

Motorola		
Phone Model	Operating System	TENS
Nexus 6	Android 7.1.1	✓
Moto X 2nd Gen	Android 6.0	✓
Moto Z	Android 7.1.1	√

One Plus		
Phone Model	Operating System	TENS
One Plus 2	Android 6.0.1	✓
One Plus 3	Android 7.0	√

Samsung		
Phone Model	Operating System	TENS
Galaxy S5	Android 6.0.1	✓
Galaxy S6	Android 7.0	✓
Galaxy S6 Edge	Android 7.0	✓
Galaxy Note 4	Android 6.0.1	✓
Galaxy Note 5	Android 7.0	✓
Galaxy S7	Android 7.0	✓
Galaxy S7 Edge	Android 7.0	✓
Galaxy S8	Android 7.0	√

Apple, iPad and iPhone are trademarks of Apple Inc., registered in the U.S. iOS is a trademark of Cisco, registered in the U.S. Android, Google, Google Play and Nexus are trademarks of Google Inc. Samsung and Galaxy S are trademarks of Samsung, registered in the U.S. LG is the registered trademark of LG Electronics, Inc. MOTOROLA is the registered trademark of Motorola Trademark Holdings, LLC. HTC and HTC One are the trademarks or registered trademarks in the U.S. and/or other countries of HTC Corporation and its affiliates.