



TENS Accessories Comparison Chart

| |  |  |  |  |  |  |
|-------------------------------|---|--|---|---|---|---|
| Model Number | Long Life Pads™ PMLLPAD | Long Life Pads™- Large PMLLPAD - L | Focus™ Band - Medium PMKBAND-M | Focus™ Band - Large PMKBAND-L | Focus™ Pad PMKNPAD | Gel Refills PMGEL |
| Reusable | Lab tested for up to 150 uses | Lab tested for up to 150 uses | Lab tested for up to 6 months | Lab tested for up to 6 months | Lab tested for up to 30 uses | Lab tested for up to 30 uses |
| Cleaning | Use tap water, air dry (up to 10 times) | Use tap water, air dry (up to 10 times) | Hand wash, air dry (up to 32 times) | Hand wash, air dry (up to 32 times) | Use moistened cloth, do not use paper cloth | Cannot be cleaned |
| Size | 4" x 2.5" each 1 Pair | 4.5" x 3" each 1 Pair | 12.5" x 15" each 1 Band | 15" x 17.7" each 1 Band | 5.1" x 2.4" each 1 Pad | 1.6" x 1.6" each 3 Pair |
| Pre-gelled, Self-adhesive | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Soft, flexible Material | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Single Person Use | Yes | Yes | Yes | Yes | Yes | Yes |
| Compatible With | PM3029, PM3030, PM3031, PM3032, PM400, PM500 | PM3029, PM3030, PM3031, PM3032, PM400, PM500 | PM710-M | PM710-L | PM710-M, PM710-L | PM311, PM800 |
| Pad/Gel Subscription Discount | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Pad Holder | ✓ | ✓ | | | | |
| Warranty | 30 day limited warranty | 30 day limited warranty | 30 day limited warranty | 30 day limited warranty | 30 day limited warranty | 30 day limited warranty |