






<p><b>Pad refills</b></p>					
<p><b>MODEL NUMBER</b></p>	<p><b>Long Life Pads™ PMLLPAD</b></p>	<p><b>Large Long Life Pads™ PMLLPAD-L</b></p>	<p><b>Gel Refills PMGEL</b></p>	<p><b>Wireless Pad - Medium PMWPAD-M</b></p>	<p><b>Wireless Pad - Large PMWPAD-L</b></p>
<p><b>Works with Omron TENS models</b></p>	<p>PM3030, PM3031, PM3032</p>	<p>PM3030, PM3031, PM3032</p>	<p>PM311</p>	<p>PM601</p>	<p>PM601</p>
<p><b>FEATURES / BENEFITS</b></p>					
<p><b>Re-usable</b></p>	<p>Lab tested for up to 150 uses</p>	<p>Lab tested for up to 150 uses</p>	<p>Lab tested for up to 30 uses</p>	<p>Lab tested for up to 30 uses</p>	<p>Lab tested for up to 30 uses</p>
<p><b>Washable</b></p>	<p>Yes, up to 10 times (use tap water, air dry)</p>	<p>Yes, up to 10 times (use tap water, air dry)</p>	<p>No</p>	<p>Yes, up to 10 times (use moistened cloth; do not use paper cloth)</p>	<p>Yes, up to 10 times (use moistened cloth; do not use paper cloth)</p>
<p><b>Soft, flexible material</b></p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p><b>Size</b></p>	<p>4" H x 2.5" W each 1 Pair</p>	<p>4.75" H x 3" W each 1 Pair</p>	<p>1" H x 1.5" W each 3 Pair</p>	<p>3.1" H x 7" W 1 pair</p>	<p>3.3" H x 8.6" W 1 pair</p>
<p><b>Pre-gelled, self- adhesive</b></p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p><b>Single person use</b></p>	<p>Yes <i>(Pad Refills available online)</i></p>	<p>Yes <i>(Pad Refills available online)</i></p>	<p>Yes <i>(Gel Refills available online)</i></p>	<p>Yes <i>(Pad Refills available online)</i></p>	<p>Yes <i>(Pad Refills available online)</i></p>
<p><b>Pad/Gel Subscription Discount</b></p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p><b>Pad Holder</b></p>		<p>✓</p>			
<p><b>Warranty</b></p>	<p>30 days limited</p>	<p>30 days limited</p>	<p>30 days limited</p>	<p>30 days limited</p>	<p>30 days limited</p>