




TENS Accessories Comparison Chart

								
Model Number	Long Life Pads™ PMLLPAD	Long Life Pads™ - Large PMLLPAD - L	Avail™ Wireless Pad - Medium PMWPAD - M	Avail™ Wireless Pad - Large PMWPAD - L	Focus™ Band - Medium PMKBAND-M	Focus™ Band - Large PMKBAND-L	Focus™ Pad PMKNPAD	Gel Refills PMGEL
Reusable	Lab tested for up to 150 uses	Lab tested for up to 150 uses	Lab tested for up to 30 uses	Lab tested for up to 30 uses	Lab tested for up to 6 months	Lab tested for up to 6 months	Lab tested for up to 30 uses	Lab tested for up to 30 uses
Cleaning	Use tap water, air dry (up to 10 times)	Use tap water, air dry (up to 10 times)	Use moistened cloth (up to 2 times); do not use paper cloth	Use moistened cloth (up to 2 times); do not use paper cloth	Hand wash, air dry (up to 32 times)	Hand wash, air dry (up to 32 times)	Use moistened cloth, do not use paper cloth	Cannot be cleaned
Size	4" x 2.5" each 1 Pair	4.5" x 3" each 1 Pair	3.1" x 7" each 1 Pad	3.3" x 8.6" each 1 Pad	12.5" x 15" each 1 Band	15" x 17.7" each 1 Band	5.1" x 2.4" each 1 Pad	1.6" x 1.6" each 3 Pair
Pre-gelled, Self-adhesive	✓	✓	✓	✓	✓	✓	✓	✓
Soft, flexible Material	✓	✓	✓	✓	✓	✓	✓	✓
Single Person Use	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Compatible With	PM3029, PM3030, PM3031, PM3032, PM400, PM500	PM3029, PM3030, PM3031, PM3032, PM400, PM500	PM601	PM601	PM710-M	PM710-L	PM710-M, PM710-L	PM311, PM800
Pad/Gel Subscription Discount	✓	✓	✓	✓	✓	✓	✓	✓
Pad Holder	✓	✓						
Warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty