

OMRON[®] MICRO *Air*[®]

Welcome to a world where you can breathe freely. We've mobilized the nebulizer so it can go wherever you go.

Break free from your respiratory restrictions with the Omron[®] MicroAir[®] Nebulizer. Spend a day on the lake or a whole afternoon at the museum without having to worry about getting home to take your next treatment. Portable, quiet relief will open up your world. Accurate medication delivery from a 6 oz. nebulizer means that any place can be the right place to take a treatment.

Actual
SIZE



Prescription
required for
purchase.

Benefits

- Portable
- Easy to use
- Handheld
- Accurate
- Tubeless
- Quiet and discreet
- Dispenses medication quickly

Features

- Compact, portable size
- Runs on two AA batteries
- 40 treatments or 4 hours of continuous use
- Universal mouthpiece adaptor
- Cutting-edge Vibrating Mesh Technology quickly dispenses medication

Share this with your doctor to
see if MicroAir is right for you.

easy
simple to
use & carry



small
a compact
solution



accurate
precise &
efficient delivery

