



LONG LIFE PADS™ REFILLS



MODEL NUMBER	Long Life Pads™ PM LLPAD	Large Long Life Pads™ PM LLPAD – L
FEATURES/BENEFITS		
Re-usable	Lab tested for up to 150 uses	Lab tested for up to 150 uses
Washable	Yes, up to 15 times (use tap water, air dry)	Yes, up to 15 times (use tap water, air dry)
Fits Omron models PM3030 & PM 3031	✓	✓
Soft, flexible material	✓	✓
Size	2.5" W x 4" H each 1 Pair	3" W x 4.75" H each 1 Pair
Pre-gelled, self-adhesive	✓	✓
Single person use	✓ <i>(Pad Refills available for additional family members)</i>	✓ <i>(Pad Refills available for additional family members)</i>
Why Long Life Pads?		
Good Value	Lasts long time, don't need to refill as often	Lasts long time, don't need to refill as often
For Multiple Body Pains	Lower back, arm, leg/foot, hip/thigh, joints	Larger-sized muscles such as lower back, hip/thigh
Pain Relief Pad Club Discounts	✓ store.omronhealthcare.com	✓ store.omronhealthcare.com
OTHER INFORMATION		
Quick Start Guide	✓	✓
Pad Holder		✓
Warranty	30 days limited	30 days limited